

General

Name of Run	231010_Q2_Run2R1.run	Name of Analysis	NOx_P04.ANL
Date of Report	10/16/2023	System	AA500
Run Start	10/10/2023 4:30:58 PM	Run Stop	5:26:52 PM
Operator Run/Rec.	CLMR / CLMR	Software Version	AACE 8.02 SP2
Comment	Recalculate from Run231010_Q2_Run2.run r_1st Batch		

Parameters

Channel	1	3	
Method	NOx	P04	Dilution
Unit	mg/L	mg/L	
Calib. Fit	Linear	Linear	
Corr.Coeff.(r)	0.9997	0.9999	
AUFS/Gain	0.42/24	0.10/103	Data
Offset	-331	4977	Field
Sensitivity	0.3338	0.0847	'#2
Sample Limit 1	-	-	
Sample Limit 2	-	-	

Results

Pk	Cup	Type	Sample ID	Value	Value	Data
0	0	B	Initial Base	-0.002	0.001	0
1	909	P	Primer	0.198	0.198	1
2	909	D	Drift	0.199	0.199	1
3	909	C	0.200	0.199	0.200	1
4	901	C	0.100	0.102	0.100	1
5	902	C	0.050	0.052	0.050	1
6	903	C	0.025	0.025	0.026	1
7	904	C	0.010	0.011	0.011	1
8	905	C	0.005	0.003	0.004	1
9	906	C	Blank	-0.001	-0.001	1
10	0	B	Baseline	-0.002	0.001	1
11	909	H1	High	0.199	0.199	1
12	905	L1	Low	0.003	0.004	1
13	905	L1	Low	0.003	0.004	1
14	906	S	Blank	-0.001	-0.001	1
15	907	K	Blank spiked	0.150	0.148	1
16	908	QC1	CRM N03	0.026 P	0.000	1
17	1	S	69	0.041	0.009	1
18	2	S	70	0.004	0.032	1
19	3	S	71	0.005	0.007	1
20	4	S	72	0.016	0.021	1
21	5	S	73	0.010	0.042	1
22	6	K	73 spiked	0.165	0.197	1
23	908	QC1	CRM N03	0.027 P	0.000	1
24	0	B	Baseline	-0.002	0.001	1
25	7	S	74	0.006	0.015	1
26	8	S	75	0.004	0.007	1
27	9	S	76	0.001	0.005	1
28	10	S	77	-0.002	0.012	1
29	11	S	78	0.001	0.014	1
30	12	S	81	0.001	0.016	1
31	13	K	81 spiked	0.028	0.164	1
32	908	QC1	CRM N03	0.029 F	0.000	1
33	0	B	Baseline	-0.002	0.001	1
34	921	D	Drift	-0.002	0.001	1
35	0	B	Final Base	-0.002	0.001	1

Corrections

Channel	1	3
Baseline	done	done
Drift	none	none
Carryover	done	done
%	0.20	0.31

Spike Recovery


Blank spiked Added	0.150	0.150
Blank spiked Found	0.151	0.148
Recovery %	100.594	98.861
73 spiked Added	0.150	0.149
73 spiked Found	0.154	0.155
Recovery %	103.111	104.004
81 spiked Added	0.150	0.150
81 spiked Found	0.027	0.148
Recovery %	18.088	99.230

Calibration

Channel 1 NOx

Type	Observed	Calculated	Target	% Error
1C	55843	0.199	0.200	-0.73
2C	30524	0.102	0.100	1.94
3C	17418	0.052	0.050	3.88
4C	10393	0.025	0.025	0.55
5C	6639	0.011	0.010	8.17
6C	4588	0.003	0.005	-40.16
7C	3446	-0.001	0.000	-


Calibration Fit: Linear
Correlation Coefficient: 0.9997
Carryover (%): 0.1995
Calibration Equation: $y = bx + a$
y = Peak Height Digit. Units
x = Concentration Unit mg/L
a = 3.803781E+003
b = 2.621142E+005



Channel 3 P04

Type	Observed	Calculated	Target	% Error
1C	60807	0.200	0.200	-0.16
2C	31951	0.100	0.100	0.27
3C	17464	0.050	0.050	0.73
4C	10322	0.026	0.025	3.05
5C	5973	0.011	0.010	7.80
6C	3906	0.004	0.005	-26.81
7C	2696	-0.001	0.000	-

Calibration Fit: Linear
Correlation Coefficient: 0.9999
Carryover (%): 0.3067
Calibration Equation: $y = bx + a$
y = Peak Height Digit. Units
x = Concentration Unit mg/L
a = 2.843657E+003
b = 2.902908E+005



Legend

* ... Sample offscale
+ ... Result higher than sample limit
- ... Result lower than sample limit
P ... Standard passed
F ... Standard failed
N ... Value not calculated or not used
R ... Resample after offscale
M ... Peak marker moved manually
d ... Pre-Diluted sample
D ... Post-Diluted sample
H ... Dual calibration used
Z ... Report negative values as zero