

General

Name of Run	240923_Q6_NP3.run	Name of Analysis	NOx_P04.ANL
Date of Report	9/23/2024	System	AA500
Run Start	9/23/2024 4:31:30 PM	Run Stop	6:52:52 PM
Operator Run/Rec. CLMR		Software Version	AACE 8.02 SP2
Comment			

Parameters

Channel	1	3	
Method	NOx	P04	Dilution
Unit	mg/L	mg/L	
Calib. Fit	Linear	Linear	
Corr.Coeff.(r)	1.0000	0.9999	
AUFS/Gain	0.24/42	0.10/99	Data
Offset	166	4893	Field
Sensitivity	0.1811	0.0851	'#2
Sample Limit 1	-	-	
Sample Limit 2	-	-	

Results

Pk	Cup	Type	Sample ID	Value	Value	Data
0	0	B	Initial Base	-0.003	-0.004	0
1	921	P	Primer	0.214	0.198	1
2	921	D	Drift	0.199	0.199	1
3	921	C	0.200	0.199	0.199	1
4	901	C	0.100	0.101	0.101	1
5	902	C	0.050	0.050	0.051	1
6	903	C	0.025	0.025	0.024	1
7	904	C	0.010	0.009	0.009	1
8	905	C	0.005	0.004	0.006	1
9	906	C	Blank	0.000	-0.001	1
10	0	B	Baseline	-0.003	-0.004	1
11	921	H1	High	0.195	0.199	1
12	905	L1	Low	0.002	0.006	1
13	905	L1	Low	0.005	0.006	1
14	906	S	Blank	-0.001	-0.001	1
15	907	K	Blank spiked	0.138	0.147	1
16	908	QC1	CRM 7/15 [0.15/0.03]	0.132	0.029	1
17	909	QC1	CRM 4/30 [0.15/0.03]	0.134	0.028	1
18	1	S	321	0.012	0.009	1
19	2	S	399	0.014	0.006	1
20	3	S	139	0.013	0.006	1
21	4	S	464	0.041	0.018	1
22	5	S	360	0.009	0.009	1
23	6	S	529	0.009	0.009	1
24	7	S	9	0.011	0.006	1
25	8	S	295	0.010	0.006	1
26	9	S	503	0.009	0.008	1
27	10	S	607	0.017	0.010	1
28	11	K	607 spiked	0.155	0.158	1
29	908	QC1	CRM 7/15 [0.15/0.03]	0.130	0.028	1
30	909	QC1	CRM 4/30 [0.15/0.03]	0.132	0.027	1
31	0	B	Baseline	-0.003	-0.004	1
32	12	S	191	0.011	0.008	1
33	13	S	347	0.014	0.006	1
34	14	S	230	0.012	0.006	1
35	15	S	204	0.011	0.006	1
36	16	S	568	0.015	0.008	1
37	17	S	126	0.009	0.005	1
38	18	S	490	0.015	0.009	1
39	19	S	425	0.009	0.005	1
40	20	S	594	0.012	0.006	1
41	21	S	217	0.010	0.009	1
42	22	K	217 spiked	0.148	0.155	1

43	908	QC1	CRM 7/15 [0.15/0.03]	0.136	0.028	1
44	909	QC1	CRM 4/30 [0.15/0.03]	0.135	0.027	1
45	0	B	Baseline	-0.003	-0.004	1
46	23	S	48	0.007	0.006	1
47	24	S	451	0.007	0.009	1
48	25	S	152	0.013	0.011	1
49	26	S	373	0.012	0.006	1
50	27	S	516	0.008	0.005	1
51	28	S	35	0.012	0.006	1
52	29	S	74	0.008	0.004	1
53	30	S	165	0.010	0.006	1
54	31	S	22	0.012	0.005	1
55	32	S	542	0.010	0.006	1
56	33	K	542 spiked	0.173	0.153	1
57	908	QC1	CRM 7/15 [0.15/0.03]	0.155	0.028	1
58	909	QC1	CRM 4/30 [0.15/0.03]	0.157	0.026	1
59	0	B	Baseline	-0.003	-0.004	1
60	34	S	477	0.016	0.005	1
61	35	S	412	0.012	0.006	1
62	36	S	243	0.016	0.006	1
63	37	S	334	0.011	0.006	1
64	38	S	308	0.015	0.006	1
65	39	S	61	0.017	0.008	1
66	40	S	282	0.009	0.009	1
67	41	S	113	0.011	0.005	1
68	42	S	100	0.007	0.005	1
69	43	S	178	0.014	0.007	1
70	44	K	178 spiked	0.173	0.153	1
71	908	QC1	CRM 7/15 [0.15/0.03]	0.148	0.028	1
72	909	QC1	CRM 4/30 [0.15/0.03]	0.146	0.026	1
73	0	B	Baseline	-0.003	-0.004	1
74	45	S	581	0.015	0.008	1
75	46	S	87	0.008	0.006	1
76	47	S	386	0.021	0.007	1
77	48	S	269	0.013	0.006	1
78	49	S	555	0.017	0.007	1
79	50	S	438	0.013	0.008	1
80	51	K	438 spiked	0.155	0.155	1
81	908	QC1	CRM 7/15 [0.15/0.03]	0.134	0.027	1
82	909	QC1	CRM 4/30 [0.15/0.03]	0.134	0.026	1
83	0	B	Baseline	-0.003	-0.004	1
84	906	S	Blank	0.000	-0.001	1
85	906	S	Blank	0.000	-0.001	1
86	906	S	Blank	-0.001	-0.001	1
87	906	S	Blank	0.000	-0.001	1
88	906	S	Blank	-0.001	-0.001	1
89	906	S	Blank	-0.001	-0.001	1
90	906	S	Blank	-0.001	-0.001	1
91	906	S	Blank	-0.001	-0.001	1
92	906	S	Blank	-0.001	-0.001	1
93	906	S	Blank	-0.001	-0.001	1
94	907	K	Blank spiked	0.136	0.146	1
95	908	QC1	CRM 7/15 [0.15/0.03]	0.144	0.029	1
96	909	QC1	CRM 4/30 [0.15/0.03]	0.142	0.028	1
97	0	B	Baseline	-0.003	-0.004	1
98	921	D	Drift	0.199	0.199	1
99	0	B	Final Base	-0.003	-0.004	1

Corrections

Channel	1	3
Baseline	done	done
Drift	done	done
Carryover	done	done
%	12.71	0.17

Spike Recovery

Blank spiked Added	0.150	0.150
Blank spiked Found	0.139	0.148
Recovery %	92.699	98.612
607 spiked Added	0.149	0.150
607 spiked Found	0.138	0.148
Recovery %	92.537	98.898
217 spiked Added	0.150	0.150


217 spiked Found	0.138	0.146
Recovery %	92.292	97.756
542 spiked Added	0.150	0.150
542 spiked Found	0.163	0.147
Recovery %	108.622	98.107
178 spiked Added	0.150	0.150
178 spiked Found	0.159	0.147
Recovery %	106.092	97.803
438 spiked Added	0.150	0.150
438 spiked Found	0.141	0.147
Recovery %	94.507	98.466
Blank spiked Added	0.150	0.150
Blank spiked Found	0.137	0.147
Recovery %	91.248	98.236

Calibration

Channel 1 NOx

Type	Observed	Calculated	Target	% Error
1C	49638	0.199	0.200	-0.31
2C	26980	0.101	0.100	0.98
3C	15347	0.050	0.050	0.92
4C	9592	0.025	0.025	1.85
5C	5877	0.009	0.010	-6.71
6C	4736	0.004	0.005	-12.53
7C	3731	0.000	0.000	-


Calibration Fit: Linear
 Correlation Coefficient: 1.0000
 Carryover (%): 12.7052
 Calibration Equation: $y = bx + a$
 y = Peak Height Digit. Units
 x = Concentration Unit mg/L
 a = 3.728984E+003
 b = 2.302516E+005



Channel 3 P04

Type	Observed	Calculated	Target	% Error
1C	59203	0.199	0.200	-0.25
2C	32150	0.101	0.100	0.83
3C	18408	0.051	0.050	1.42
4C	11175	0.024	0.025	-2.69
5C	7018	0.009	0.010	-8.33
6C	6251	0.006	0.005	27.39
7C	4257	-0.001	0.000	-

Calibration Fit: Linear
 Correlation Coefficient: 0.9999
 Carryover (%): 0.1746
 Calibration Equation: $y = bx + a$
 y = Peak Height Digit. Units
 x = Concentration Unit mg/L
 a = 4.504625E+003
 b = 2.741785E+005



Legend

* ... Sample offscale
+ ... Result higher than sample limit
- ... Result lower than sample limit
P ... Standard passed
F ... Standard failed
N ... Value not calculated or not used
R ... Resample after offscale
M ... Peak marker moved manually
d ... Pre-Diluted sample
D ... Post-Diluted sample
H ... Dual calibration used
Z ... Report negative values as zero